

GFP-Tag Mouse Monoclonal Antibody(1G6)

Catalog TDY009C TDY009F

Tel: 010-82908854

Quantity 50 μ L 100 μ L

Free: 400-0620-621

Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	N/A	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:5,000

Optimal dilutions should be determined by the end user.

Specificity: The GFP tag antibody can recognize GFP and GFP tag fusion proteins.

Background: The **green fluorescent protein (GFP)** is a protein composed of 238 amino acid residues (26.9kDa) that exhibits bright green fluorescence when exposed to light in the blue to ultraviolet range. In cell and molecular biology, the GFP gene is frequently used as a reporter of expression.



2 μ g GFP fusion protein+ Primary antibody dilution at

1、 1:5,000

2、 1:10,000