

RFP-Tag Mouse Monoclonal Antibody(5G4)

Catalog	TDY011C	TDY011F	Tel: 010-82908854
			Free: 400-0620-621
Quantity	50μL	100μL	Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	N/A	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:5,000

Optimal dilutions should be determined by the end user.

Specificity: The RFP tag antibody can recognize RFP tag fusion proteins.

Background: The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.



1μg RFP fusion protein+ Primary antibody dilution at

1、 1:3,000

2、 1:5,000

3、 1:10,000