

KT3-Tag Mouse Monoclonal Antibody(3D11)

Catalog	TDY015C	TDY015F	Tel: 010-82908854
			Free: 400-0620-621
Quantity	50μL	100μL	Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	N/A	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:5,000

Optimal dilutions should be determined by the end user.

Specificity: The KT3 tag antibody can recognize KT3 tag fusion proteins.

Background: KT3-tag peptide KPPTPPPEPET conjugated to KLH. KT3-Tag antibody can recognize C-terminal, internal, and N-terminal KT3-tagged proteins.



1μg KT3 fusion protein+ Primary antibody dilution at

1、 1:5,000

2、 1:10,000