

CK19 Mouse Monoclonal Antibody(11F5)

Catalog TDY074C TDY074F

Tel: 010-82908854

Quantity 50μL 100μL

Free: 400-0620-621

Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,M,R	~44KD	IgG1


Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: IHC: 1:200

Optimal dilutions should be determined by the end user.

Specificity: The CK19 Mouse Monoclonal antibody detects endogenous CK19 proteins.

Background: This protein is a member of the keratin family. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the periderm, the transiently superficial layer that envelopes the developing epidermis. Keratin 19 is not expressed in hepatocytes, therefore, antibody to keratin 19 is useful in the identification of liver metastasis. The degree of keratin 19 positivity in breast cancer distinguishes malignant from benign tumours. Keratin 19 is often co-expressed with keratin 7.



IHC staining of human colon cancer tissue with CK19 mouse mAb(11F5)
diluted at 1:200.