

## ECFP Mouse Polyclonal Antibody

Catalog TDY058C TDY058F

Tel: 010-82908854

Quantity 50 $\mu$ L 100 $\mu$ L

Free: 400-0620-621

Web: www.tdybio.com

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IHC, IF	N/A	N/A	IgG

**Storage Buffer & Condition:** PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


**Recommended dilutions:** WB: 1:2,000-5,000 IHC: 1:200 IF: 1:200

**Optimal dilutions should be determined by the end user.**

**Specificity:** The ECFP Antibody can detects ECFP & CFP tag protein.

**Alternative Names:** Cyan fluorescent protein, Enhanced Cyan fluorescent protein

**Background:** Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.



Western blot analysis of ECFP Recombinant protein with TDY058 diluted at

1) 1:5,000

2) 1:10,000