

S1-Tag Mouse Monoclonal Antibody(7B10)

Catalog	TDY449C	TDY449F	Tel: 010-80117836
Quantity	50 μ L	100 μ L	Web: www.tdybio.com
			Entrez-Gene ID# , Swiss-Prot Acc.#

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	ALL	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions: WB: 1:5,000-10,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects Recombinant S1-Tag protein.

Alternative Names: N/A

Background: S1 tag is an epitope tag composed of a nine-residue peptide, NANNPDWDF, derived from the hepatitis B virus preS1 region.



Western blot analysis of Recombinant S1-Tag Protein with S1-Tag Mouse
mAb TDY449 diluted at 1:5,000.