

Endothelin B Receptor Rabbit Polyclonal Antibody(A221)

Catalog TDY522C TDY522F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# 1910, Swiss-Prot Acc.#P24530

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~50KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.


Specificity: Antibody can detects endogenous Endothelin B Receptor protein.

Alternative Names: ABCDS,Ednra,EDNRB,ET-B,ET-BR,HSCR2,WS4A

Background:



Immunohistochemical analysis of paraffin-embedded Human Prostate Tissue using Endothelin B Receptor (TDY522) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat Lung Tissue using Endothelin B Receptor (TDY522) Rabbit pAb diluted at 1:200.