



KV1.3 Rabbit Polyclonal Antibody(A256)

Catalog TDY548C TDY548F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50 μ L 100 μ L

Entrez-Gene ID# 3738, Swiss-Prot Acc.#P22001

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H,R,M	~64KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect endogenous KV1.3 protein.

Alternative Names: HGK5, KCNA3, MK3, Potassium channel 3, voltage gated potassium channel protein kv1.3



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue with KV1.3 Rabbit pAb TDY548 diluted at 1:2,000.