

## c-Jun Rabbit Polyclonal Antibody

Catalog	TDY781C	TDY781F	Tel: 010-82908854
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID#3725, Swiss-Prot Acc.#P05427

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H, M, R	43, 48KD	IgG

**Storage Buffer & Condition:**Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C. Do not aliquot the antibody.**


**Recommended dilutions:**WB: 1:1,000-2,000

**Optimal dilutions should be determined by the end user.**

**Specificity:**Antibody can detects endogenous c-Jun proteins.

**Alternative Names:**Activator Protein 1 antibody, AP1 antibody, p39 antibody

**Background:**c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1).



Western blot analysis of 1) 293T Cell Lysate, 2) Rat Kidney Tissue Lysate, 3) 3T3 Cell

using c-Jun (TDY781) Rabbit pAb diluted at 1:2000.